
^{13}C Atemtestprotokolle - Darm

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■ Empfohlene Literatur

1. Romagnuolo J, Schiller D, Bailey RJ. Using breath tests wisely in a gastroenterology practice: an evidence-based review of indications and pitfalls in interpretation. *Am J Gastroenterol.* 2002 May;97(5):1113-26. Review.
2. Vonk RJ, Stellaard F, Priebe MG, Koetse HA, Hagedoorn RE, De Bruijn S, Elzinga H, Lenoir-Wijnkoop I, Antoine JM. The ¹³C/2H-glucose test for determination of small intestinal lactase activity. *Eur J Clin Invest.* 2001 Mar;31(3):226-33.
3. Stellaard F, Koetse HA, Elzinga H, Boverhof R, Tjoonk R, Klimp A, Vegter D, Liesker J, Vonk RJ. ¹³C-carbohydrate breath tests: impact of physical activity on the rate-limiting step in lactose utilization. *Scand J Gastroenterol.* 2000 Aug;35(8):819-23.
4. Koetse HA, Vonk RJ, Pasterkamp S, Pal J, de Bruijn S, Stellaard F. Variations in colonic H₂ and CO₂ production as a cause of inadequate diagnosis of carbohydrate maldigestion in breath tests. *Scand J Gastroenterol.* 2000 Jun;35(6):607-11.
5. Vonk RJ, Lin Y, Koetse HA, Huang C, Zeng G, Elzinga H, Antoine J, Stellaard F. Lactose (mal)digestion evaluated by the ¹³C-lactose digestion test. *Eur J Clin Invest.* 2000 Feb;30(2):140-6.
6. Koetse HA, Stellaard F, Bijleveld CM, Elzinga H, Boverhof R, van der Meer R, Vonk RJ, Sauer PJ. Non-invasive detection of low-intestinal lactase activity in children by use of a combined ¹³CO₂/H₂ breath test. *Scand J Gastroenterol.* 1999 Jan;34(1):35-40.
7. Murray RD, Boutton TW, Klein PD, Gilbert M, Paule CL, MacLean WC Jr. Comparative absorption of [¹³C]glucose and [¹³C]lactose by premature infants. *Am J Clin Nutr.* 1990 Jan;51(1):59-66.
8. Hiele M, Ghoois Y, Rutgeerts P, Vantrappen G, Carchon H, Eggermont E. ¹³CO₂ breath test using naturally ¹³C-enriched lactose for detection of lactase deficiency in patients with gastrointestinal symptoms. *J Lab Clin Med.* 1988 Aug;112(2):193-200.
9. MacLean WC Jr, Fink BB, Schoeller DA, Wong W, Klein PD. Lactose assimilation by full-term infants: relation of [¹³C] and H₂ breath tests with fecal [¹³C] excretion. *Pediatr Res.* 1983 Aug;17(8):629-33.

2. Orocoecale Transitzeit - ¹³C-Laktose Ureid Atemtest

■ Empfohlene Literatur

1. Christian M, Morrison D, Dodson B, Preston T, Amarri S, Franchini F, Edwards C, Weaver L. Measurement of oro-cecal transit time in young children using lactose [¹³C] ureide requires further validation. *J Pediatr Gastroenterol Nutr.* 2002 May;34(5):570-1; author reply 571.
2. Van Den Driessche M, Van Malderen N, Geypens B, Ghooos Y, Veereman-Wauters G. Lactose-[¹³C]ureide breath test: a new, noninvasive technique to determine orocecal transit time in children. *J Pediatr Gastroenterol Nutr.* 2000 Oct;31(4):433-8.
3. Jackson AA, Bundy R, Hounslow A, Murphy JL, Wootton SA. Metabolism of lactose-[¹³C]ureide and lactose-[¹⁵N, ¹⁵N]ureide in normal adults consuming a diet marginally adequate in protein. *Clin Sci (Lond).* 1999 Nov;97(5):547-55.
4. Geypens B, Bennink R, Peeters M, Evenepoel P, Mortelmans L, Maes B, Ghooos Y, Rutgeerts P. Validation of the lactose-[¹³C]ureide breath test for determination of orocecal transit time by scintigraphy. *J Nucl Med.* 1999 Sep;40(9):1451-5.
5. Morrison DJ, Dodson B, Preston T. Measurement of urinary total ¹³C and ¹³C urea by isotope ratio mass spectrometry after administration of lactose [¹³C]-ureide. *Rapid Commun Mass Spectrom.* 1999;13(13):1252-6.
6. Morrison DJ, Zavoshy R, Edwards CA, Dodson B, Preston T, Weaver LT. Lactose [¹³C]Ureide as a marker for colonic fermentation and the deconvolution of a complex ¹³CO₂ breath test curve. *Biochem Soc Trans.* 1998 May;26(2):S184.
7. Ruemmele FM, Heine WE, Keller KM, Lentze MJ. Metabolism of glycosyl ureides by human intestinal brush border enzymes. *Biochim Biophys Acta.* 1997 Aug 29;1336(2):275-80.
8. Wutzke KD, Heine WE, Plath C, Leitzmann P, Radke M, Mohr C, Richter I, Gulzow HU, Hobusch D. Evaluation of oro-coecal transit time: a comparison of the lactose-[¹³C, ¹⁵N]ureide ¹³CO₂- and the lactulose H₂-breath test in humans. *Eur J Clin Nutr.* 1997 Jan;51(1):11-9.
9. Heine WE, Berthold HK, Klein PD. A novel stable isotope breath test: ¹³C-labeled glycosyl ureides used as noninvasive markers of intestinal transit time. *Am J Gastroenterol.* 1995 Jan;90(1):93-8.

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■ Empfohlene Literatur

1. Dellert SF, Nowicki MJ, Farrell MK, Delente J, Heubi JE. The ¹³C-xylose breath test for the diagnosis of small bowel bacterial overgrowth in children. *J Pediatr Gastroenterol Nutr.* 1997 Aug;25(2):153-8.
2. Lim H, Wagner DA, Toskes PP. ¹³C-xylose breath test for bacterial overgrowth. *Gastroenterology* 1993 April;104(4): A259
3. Valdovinos MA, Camilleri M, Thomforde GM, Frie C. Reduced accuracy of ¹⁴C-D-xylose breath test for detecting bacterial overgrowth in gastrointestinal motility disorders. *Scand J Gastroenterol.* 1993 Nov;28(11):963-8.
4. Chang CS, Chen GH, Kao CH, Wang SJ, Peng SN, Huang CK, Poon SK. Increased accuracy of the carbon-¹⁴ D-xylose breath test in detecting small-intestinal bacterial overgrowth by correction with the gastric emptying rate. *Eur J Nucl Med.* 1995 Oct;22(10):1118-22.
5. Stotzer PO, Kilander AF. Comparison of the 1-gram (¹⁴C)-D-xylose breath test and the 50-gram hydrogen glucose breath test for diagnosis of small intestinal bacterial overgrowth. *Digestion.* 2000;61(3):165-71.
6. King CE, Toskes PP, Guilarte TR, Lorenz E, Welkos SL. Comparison of the one-gram d-[¹⁴C]xylose breath test to the [¹⁴C]bile acid breath test in patients with small-intestine bacterial overgrowth. *Dig Dis Sci.* 1980 Jan;25(1):53-8.
7. King CE, Toskes PP. Breath tests in the diagnosis of small intestine bacterial overgrowth. *Crit Rev Clin Lab Sci.* 1984;21(3):269-81. Review.
8. Lewis SJ, Young G, Mann M, Franco S, O'Keefe SJ. Improvement in specificity of [¹⁴C]d-xylose breath test for bacterial overgrowth. *Dig Dis Sci.* 1997 Aug;42(8):1587-92.
9. Rumessen JJ. [¹⁴C]D-xylose breath test for small intestinal bacterial overgrowth. *Gastroenterology.* 1989 Jan;96(1):273-4.
10. Riordan SM, Mclver CJ, Duncombe VM, Bolin TD, Thomas MC. Factors influencing the 1-g ¹⁴C-D-xylose breath test for bacterial overgrowth. *Am J Gastroenterol.* 1995 Sep;90(9):1455-60.
11. King CE, Toskes PP. Comparison of the 1-gram [¹⁴C]xylose, 10-gram lactulose-H₂, and 80-gram glucose-H₂ breath tests in patients with small intestine bacterial overgrowth. *Gastroenterology.* 1986 Dec;91(6):1447-51.
12. Schneider A, Novis B, Chen V, Leichtman G. Value of the ¹⁴C-D-xylose breath test in patients with intestinal bacterial overgrowth. *Digestion.* 1985;32(2):86-91.
13. King CE, Toskes PP, Spivey JC, Lorenz E, Welkos S. Detection of small intestine bacterial overgrowth by means of a ¹⁴C-D-xylose breath test. *Gastroenterology.* 1979 Jul;77(1):75-82.

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© Campro Scientific GmbH