

Protocols - ^{13}C Breath Tests - LIVER




Table of Contents

Introduction.....	3
• ¹³ C Breath Test protocols	3
• ¹³ C Breath Testing: principle and requirements	3
• Preparation of the patient.....	3
• Administration of ¹³ C labeled substrate	3
• Collection of breath samples.....	3
• Measurement of ¹³ C enrichment	3
• Calculation of the end result.....	3
• Applications.....	4
• Literature	4
• Note.....	4
• Authors.....	4
Test Protocols Liver.....	5
• Microsomal Liver Function	5
¹³ C-Aminopyrine Breath Test	
• Microsomal Liver Function	6
¹³ C-Methacetin Breath Test	
• Mitochondrial Liver Funktion	7
¹³ C-Ketoisocaproic Acid Breath Test (¹³ C-KICA)	
• Cytosolic Liver Funktion.....	8
¹³ C-Galactose Breath Test	
Literature	9
Liver	9
• Microsomal Liver Function	9
¹³ C-Aminopyrine Breath Test	
• Microsomal Liver Function	10
¹³ C-Methacetin Breath Test	
• Mitochondrial Liver Funktion	11
¹³ C-Ketoisocaproic Acid Breath Test (¹³ C-KICA)	
• Cytosolic Liver Funktion	12
¹³ C-Galactose Breath Test	

Introduction

■ **¹³C Breath Test protocols**

This folder contains a set of protocols describing the principles and general test procedures for today's most relevant ¹³C Breath Tests to study specific functions of the liver, pancreas, stomach and intestine. The list will be updated regularly adding additional tests or additional information on the already described tests.

The information is meant as a start to enter the field of stable isotope ¹³C Breath testing initiated by the interest in a specific test. To actually be able to introduce a test in your hospital you must familiarize yourself with basic knowledge of breath testing with ¹³C substrates and the existing knowledge on the particular application of interest. There is no such thing as a standard protocol for all tests.

■ **¹³C Breath Testing: principle and requirements**

A ¹³C Breath Test consists of the administration to a patient of a ¹³C labeled substrate that is metabolized by a specific enzyme system resulting in ¹³CO₂ as the end product. To monitor the enzyme response ¹³C enrichment is measured in breath CO₂.

The total procedure of ¹³C Breath testing includes the definition of the preparation of the patient before the test, administration of the ¹³C labeled substrate, collection of breath samples, measurement of ¹³C enrichment in breath CO₂ and the calculation of an end result.

■ **Preparation of the patient**

In general, tests will be performed in the fasting state and the patient should be at a low and stable level of natural ¹³C abundance. Therefore, the patient must be instructed to avoid eating ¹³C enriched foods such as corn products, cane sugar, pineapple and tequila the last days before the test and to come to the clinic fasted. In certain cases (¹³C Lactose-Ureide breath test) the patient must be pretreated with unlabeled substrate to stimulate the involved enzyme system.

■ **Administration of ¹³C labeled substrate**

The test substrate may be administered as a simple solution in water with or without a standardized test meal. Sometimes it needs to be incorporated into a specific ingredient of the meal. The test meal and the dose of substrate may be different for adults and children.

■ **Collection of breath samples**

Every protocol has its own time schedule of breath collections. The number of samples may be as small as 2 or more than 20. To define the ¹³C enrichment in breath CO₂, it is also necessary to obtain at least two breath samples before the ingestion of the ¹³C substrate to determine the natural background of ¹³C abundance. The methodology of collecting breath samples is dependent on the technology to determine the ¹³C enrichment. The protocols are based on Continuous Flow Isotope Ratio Mass Spectrometry as the analytical technique. In this case breath samples are simply blown through a straw into special 10 ml gas collection tubes that directly fit into the sample tray of the instruments. In case of Infrared technology special bags provided by the instrument manufacturer must be used.

■ **Measurement of ¹³C enrichment**

To determine the ¹³C abundance in breath CO₂ you need the availability of Isotope Ratio Mass Spectrometry (IRMS) or specialized Infrared instrumentation. The protocols are based on Isotope Ratio Mass Spectrometry. For a number of tests (Aminopyrine, Methacetin, Urea) Infrared Spectroscopy has proven to be a valid alternative analytical technique. For other tests Infrared technology has not yet been validated so far. In principle the test substrate is not a determinant of the validity of the analytical technique. It is the level of ¹³C enrichment that determines the analytical requirement. Validation of Infrared analysis for other application is recommended, as it is recommended to validate any breath test in your own clinical laboratory. You may have instrumentation available or contact a service center for the analyses.

■ **Calculation of the end result**

For some tests the only calculation needed is the subtraction of the natural background value from the measured value at a defined time. In other cases it is necessary to calculate the amount of ¹³C that is recovered in breath during the experimental period. In a third type of application the time course of the enrichment appearance is of importance requiring calculation of the rate of appearance.

■ Applications

In the present update the following tests have been described:

Table of contents of applications

Function		¹³ C Substrate	Page
Liver			
1.	Microsomal (cirrhosis)	¹³ C-Aminopyrine	5
2.	Microsomal	¹³ C-Methacetin	6
3.	Mitochondrial	¹³ C-Ketoisocaproic Acid	7
4.	Cytosolic	¹³ C-Galactose	8
Literature			9

■ Literature

Included is a list of literature references that will introduce you to the most important articles describing aspects of the different tests described in the protocols.

■ Note

Great care has been taken over the composition of the text, figures and tables. The possibility of errors however, cannot be excluded completely. Therefore Campro Scientific GmbH and the authors cannot accept any legal or other liability with respect to incorrect details and their consequences. The authors would be grateful to receive suggestions for improvements and information about errors.

This information may not be photocopied, duplicated or translated into another language neither partially nor completely, in any form without the written agreement of Campro Scientific GmbH.

■ Authors

1. Dr. F. Stellaard
University Hospital Groningen, the Netherlands
Dept. Pathology and Laboratory Medicine
And Centre for Liver, Intestinal and Metabolic Disease

2. Dr. Ahmad Rajabi
Campro Scientific GmbH
Berlin, Germany

© Campro Scientific GmbH

European Headquarters

P.O. Box 37 03 31
D-14133 Berlin
Germany
Tel. : +49(0)30.629.01.89.0
Fax : +49(0)30.629.01.89.89
E-mail : info@campro.eu
Web : www.campro.eu

Dutch Sales Office

P.O. Box 316
NL-3900 AH Veenendaal
The Netherlands
Tel. : +31(0)318.529.437
Fax : +31(0)318.542.181
E-mail : info@campro.eu
Web : www.campro.eu

Please contact us for technical and price information.

Test Protocols Liver

1. Microsomal Liver Function ¹³C-Aminopyrine Breath Test

■ Principle

¹³C₂-Aminopyrine contains two methyl groups both labeled with the non-radioactive isotope ¹³C. After oral administration, ¹³C₂-Aminopyrine is completely absorbed and cleared by the liver.

¹³C₂-Aminopyrine is demethylated by the microsomal cytochrome P450 dependent N-demethylase enzyme. By this reaction ¹³C-Formaldehyde is formed of which 50% is oxidized to bicarbonate. The degree of appearance of ¹³C in breath CO₂ reflects the degree of demethylation.

■ Applicability of ¹³C-Aminopyrine Breath Test

¹³C-Aminopyrine Breath Test has so far been applied to adults, children and even neonates.

■ Applications

¹³C-Aminopyrine Breath Test is used to detect diminished microsomal liver function by comparison with a control range. The result reflects microsomal mass and is independent of liver blood flow.

The ¹³C-Aminopyrine Breath Test, most useful in cirrhotic patients, is correlated with severity of liver disease and promoted as a reliable method for the optimal timing for liver transplantation. Also P450 dependent enzyme induction by therapeutic means can be monitored by ¹³C-Aminopyrine Breath Test.

■ Protocol

Adults: The ¹³C-Aminopyrine Breath Test is performed after an overnight fast. A dose of 2 mg/kg body weight ¹³C₂-Aminopyrine is administered orally after dissolution in about 100 ml water. Breath samples are collected before (2x) and every 30 minutes for 120 minutes (2 h) after the ingestion of the ¹³C₂-Aminopyrine. ¹³C enrichment in breath CO₂ is determined by Isotope Ratio Mass Spectrometry (IRMS). The cumulative percentage of ¹³C recovered in breath during the 120 minutes, collection period is used as the diagnostic parameter.

Children: A dose of 2-5 mg/kg body weight is used for children. The same time schedule for breath collections can be used.

■ Interpretation of test results

It is advised to obtain your own internal control values. Generally a cut-off value of 7% for the 120 minutes cumulative recovery may serve as a starting point in the case of adults and older children. Values below 7% are indicate diminished microsomal function. Within the first year of life the microsomal mass develops and the test results are therefore age dependent.

■ Precautions

Microsomal cytochrome P450 enzyme activity may be induced by alcohol or drugs (spironolactone, glutethimide, diphenylhydantoin) or depressed (cimetidine) which affects the outcome of the test. In exceptional cases chronic dosing of aminopyrine has been associated with aminopyrine-induced agranulocytosis.

■ Summary

	Dose	Samples	
Adults	2 mg/kg body weight ¹³ C ₂ -Aminopyrine	2	Before administration
		4	Every 30 minutes for 120 minutes after administration (2 hours)
Children	2-5 mg/kg body weight ¹³ C ₂ -Aminopyrine	2	Before administration
		4	Every 30 minutes for 120 minutes after administration (2 hours)

2. Microsomal Liver Function

¹³C-Methacetin Breath Test

■ Principle

¹³C-Methacetin contains one methyl group labeled with the non-radioactive isotope ¹³C. After oral administration ¹³C-Methacetin is completely absorbed and cleared by the liver. ¹³C-Methacetin is demethylated by the microsomal cytochrome P450 dependent N-demethylase enzyme. By this reaction ¹³C-formaldehyde is formed of which 50% is oxidized to bicarbonate. The degree of appearance of ¹³C in breath CO₂ reflects the degree of demethylation.

■ Applicability of ¹³C-Methacetin Breath Test

¹³C-Methacetin Breath Test has so far been applied to adults and children.

■ Applications

¹³C-Methacetin Breath Test is used to detect diminished microsomal liver function by comparison with a control range. The result reflects microsomal mass and is independent of liver blood flow. The ¹³C-Methacetin Breath Test, most useful in cirrhotic patients, is correlated with severity of liver disease and promoted as a reliable method for the optimal timing for liver transplantation. Also P450 dependent enzyme induction by therapeutic means can be monitored by the ¹³C-Methacetin Breath Test.

■ Protocol

Adults: The ¹³C-Methacetin Breath Test is performed after an overnight fast. A dose of 2 mg/kg body weight ¹³C-Methacetin is administered orally after dissolution in about 100 ml water. Breath samples are collected before (2x) and every 10 minutes for 30 minutes after ingestion of the ¹³C-Methacetin. ¹³C enrichment in breath CO₂ is determined by Isotope Ratio Mass Spectrometry (IRMS). The cumulative percentage of ¹³C recovered in breath during the 30 minutes collection period is used as the diagnostic parameter.

■ Interpretation of test results

It is advised to obtain your own internal control values. Generally a cut-off value of 8% for the 30 minutes cumulative recovery may serve as a starting point.

■ Precautions

Microsomal P450 enzyme activity may be induced by alcohol or drugs (spironolactone, glutethimide, diphenylhydantoin) or depressed (cimetidine), which affects the outcome of the test.

■ Summary

	Dose	Samples	
	2 mg/kg body weight ¹³ C-Methacetin	2	Before administration
		3	Every 10 minutes for 30 minutes after administration (0.5 h)

3. Mitochondrial Liver Funktion

¹³C-Ketoisocaproic Acid Breath Test (¹³C-KICA)

■ Principle

Ketoisocaproic Acid (KICA) is the result of branched chain amino acid aminotransferase degradation of leucine. ¹³C-Ketoisocaproic acid (¹³C-KICA) is labeled with the non-radioactive isotope ¹³C. After oral administration, ¹³C-KICA is completely absorbed and cleared by the liver. ¹³C-KICA is decarboxylated almost exclusively in the hepatic mitochondria producing ¹³CO₂. The degree of appearance of ¹³C in breath CO₂ reflects the degree of mitochondrial decarboxylation.

■ Applicability of ¹³C-KICA Breath Test

¹³C-KICA Breath Test has so far been applied to adults.

■ Applications

¹³C-KICA Breath Test is used to detect diminished mitochondrial liver function by comparison with a control range. The result reflects not a general but a specific ethanol-related impairment of mitochondrial function and is used to determine the involvement of ethanol abuse in liver disease.

■ Protocol

The ¹³C-KICA Breath Test is performed after an overnight fast. A dose of 1 mg/kg body weight ¹³C KICA is administered orally together with 20 mg/kg L-Leucine after dissolution in 200 ml 0.1N citric acid solution or 200 ml orange juice. The L-Leucine is added to increase the decarboxylation of KICA. Breath samples are collected before (2x) and 10, 15, 20, 25, 30 and 60 minutes after administration of the ¹³C-KICA. ¹³C enrichment in breath CO₂ is determined by Isotope Ratio Mass Spectrometry (IRMS). The cumulative percentage of ¹³C recovered in breath during the 60 minutes collection period is used as the diagnostic parameter.

■ Interpretation of test results

It is advised to obtain your own internal control values. Generally a cut-off value of 7% for the 60 minutes cumulative recovery may serve as a starting point. Values below 7% indicate for impaired hepatic mitochondrial function.

■ Precautions

No contraindications or side effects have been described so far.

■ Summary

	Dose	Samples	
	1 mg/kg body weight ¹³ C -KICA + 20 mg/kg body weight L-Leucine	2	Before administration
		6	10, 15, 20, 25, 30 and 60 minutes after administration (1 h)

4. Cytosolic Liver Funktion ¹³C-Galactose Breath Test

■ Principle

Galactose is metabolized primarily in the liver resulting in glucose. The rate limiting step is the action of galactose kinase. After oral administration, (1-¹³C)-Galactose is completely absorbed and cleared by the liver. ¹³C-Galactose is converted almost exclusively in the hepatic cytosol producing ¹³C-Glucose. To a large extent ¹³C-Glucose is oxidized to ¹³CO₂. The degree of appearance of ¹³C in breath CO₂ reflects the degree of cytosolic conversion of galactose to galactose-1-phosphate.

■ Applicability of ¹³C-Galactose Breath Test

¹³C-Galactose Breath Test has so far been applied to adults. Preliminary experiments have been performed in children.

■ Applications

¹³C-Galactose Breath Test is used to detect diminished cytosolic liver function by comparison with a control range. The result reflects the severity of liver fibrosis in chronic hepatitis C and is proposed as a prognostic factor in the follow-up of chronic hepatitis C.

■ Protocol

The ¹³C-Galactose Breath Test is performed after an overnight fast. A dose of 5 mg/kg body weight (1-¹³C)-Galactose is administered orally together with 495 mg/kg unlabeled galactose as a 25% solution in water. Breath samples are collected before (2x) and every 10 minutes for 1 hour after the administration of the ¹³C-Galactose. ¹³C enrichment in breath CO₂ is determined by Isotope Ratio Mass Spectrometry (IRMS). The cumulative percentage of ¹³C recovered in breath / hour during the 60 minutes collection period after ingestion of (1-¹³C)-Galactose, is used as the diagnostic parameter.

■ Interpretation of test results

It is advised to obtain own internal control values. Generally a cut-off value of 4% for the 60 minutes recovery of ¹³C / hour may serve as a starting point. Values below 4% indicate the presence of fibrosis.

■ Precautions

Alcohol consumption, galactosemia and diabetes may disturb the outcome of the ¹³C-Galactose Breath Test.

■ Summary

	Dose	Samples	
Adults	5 mg/kg body weight (1- ¹³ C)-Galactose + 495 mg/kg unlabeled Galactose	2	Before administration
		6	Every 10 minutes for 60 minutes after administration (1 hour)

Literature

Liver

1. Microsomal Liver Function ¹³C-Aminopyrine Breath Test

■ Recommended literature

1. Herold C, Ganslmayer M, Ocker M, Zopf S, Gailer B, Hahn EG, Schuppan D. Inducibility of microsomal liver function may differentiate cirrhotic patients with maintained compared with severely compromised liver reserve. *J Gastroenterol Hepatol*. 2003 Apr;18(4):445-449.
2. Herold C, Regn S, Ganslmayer M, Ocker M, Hahn EG, Schuppan D. Can quantitative tests of liver function discriminate between different etiologies of liver cirrhosis? *Dig Dis Sci*. 2002 Dec;47(12):2669-73.
3. Caubet MS, Laplante A, Caille J, Brazier JL. [¹³C]aminopyrine and [¹³C]caffeine breath test: influence of gender, cigarette smoking and oral contraceptives intake. *Isotopes Environ Health Stud*. 2002 Jun;38(2):71-7.
4. Herold C, Ganslmayer M, Deynet C, Hahn EG, Schuppan D. Quantitative testing of liver function compared to prognostic scores in patients with primary biliary cirrhosis. *Liver*. 2002 Apr;22(2):159-65.
5. Giannini E, Fasoli A, Chiarbonello B, Malfatti F, Romagnoli P, Botta F, Testa E, Polegato S, Fumagalli A, Testa R. ¹³C-aminopyrine breath test to evaluate severity of disease in patients with chronic hepatitis C virus infection. *Aliment Pharmacol Ther*. 2002 Apr;16(4):717-25.
6. Fasoli A, Giannini E, Botta F, Romagnoli P, Risso D, Celle G, Testa R. ¹³CO₂ excretion in breath of normal subjects and cirrhotic patients after ¹³C-aminopyrine oral load. Comparison with MEGX test in functional differentiation between chronic hepatitis and liver cirrhosis. *Hepatogastroenterology*. 2000 Jan-Feb;47(31):234-8.
7. Opekun AR Jr, Klein PD, Graham DY. [¹³C]Aminopyrine breath test detects altered liver metabolism caused by low-dose oral contraceptives. *Dig Dis Sci*. 1995 Nov;40(11):2417-22.
8. Urbain D, Muls V, Thys O, Ham HR. Aminopyrine breath test improves long-term prognostic evaluation in patients with alcoholic cirrhosis in Child classes A and B. *J Hepatol*. 1995 Feb;22(2):179-83.
9. Mion F, Queneau PE, Rousseau M, Brazier JL, Paliard P, Minaire Y. Aminopyrine breath test: development of a ¹³C-breath test for quantitative assessment of liver function in humans. *Hepatogastroenterology*. 1995 Nov-Dec;42(6):931-8.
10. Villeneuve JP, Infante-Rivard C, Ampelas M, Pomier-Layrargues G, Huet PM, Marleau D. Prognostic value of the aminopyrine breath test in cirrhotic patients. *Hepatology*. 1986 Sep-Oct;6(5):928-31.
11. Shulman RJ, Irving CS, Boutton TW, Wong WW, Nichols BL, Klein PD. Effect of infant age on aminopyrine breath test results. *Pediatr Res*. 1985 May;19(5):441-5.
12. Jager-Roman E, Rating D, Platzek T, Helge H. Development of N-demethylase activity measured with the ¹³C-aminopyrine breath test. *Eur J Pediatr*. 1982 Oct;139(2):129-34

2. Microsomal Liver Function

¹³C-Methacetin Breath Test

■ Recommended literature

1. Ciccocioppo R, Candelli M, Di Francesco D, Ciocca F, Taglieri G, Armuzzi A, Gasbarrini G, Gasbarrini A. Study of liver function in healthy elderly subjects using the ¹³C-methacetin breath test. *Aliment Pharmacol Ther.* 2003 Jan;17(2):271-7.
2. Lara Baruque S, Razquin M, Jimenez I, Vazquez A, Gisbert JP, Pajares JM. ¹³C-phenylalanine and ¹³C-methacetin breath test to evaluate functional capacity of hepatocyte in chronic liver disease. *Dig Liver Dis.* 2000 Apr;32(3):226-32.
3. Adamek RJ, Goetze O, Boedeker C, Pfaffenbach B, Luybaerts A, Geypens B. ¹³C-methacetin breath test: isotope-selective nondispersive infrared spectrometry in comparison to isotope ratio mass spectrometry in volunteers and patients with liver cirrhosis. *Z Gastroenterol.* 1999 Dec;37(12):1139-43.
4. Pfaffenbach B, Gotze O, Szymanski C, Hagemann D, Adamek RJ. [The ¹³C-methacetin breath test for quantitative noninvasive liver function analysis with an isotope-specific nondispersive infrared spectrometer in liver cirrhosis]. *Dtsch Med Wochenschr.* 1998 Dec 4;123(49):1467-71.
5. Klatt S, Taut C, Mayer D, Adler G, Beckh K. Evaluation of the ¹³C-methacetin breath test for quantitative liver function testing. *Z Gastroenterol.* 1997 Aug;35(8):609-14.
6. Iikura Y, Iwasaki A, Tsubaki T, Akasawa A, Onda T, Katsunuma T, Miura K, Ebisawa M, Saito H, Koya N, et al. Study of liver function in infants with atopic dermatitis using the ¹³C-methacetin breath test. *Int Arch Allergy Immunol.* 1995 May-Jun;107(1-3):189-93.
7. Matsumoto K, Suehiro M, Iio M, Kawabe T, Shiratori Y, Okano K, Sugimoto T. [¹³C]methacetin breath test for evaluation of liver damage. *Dig Dis Sci.* 1987 Apr;32(4):344-8
8. Powers HJ. The potential of a ¹³C-methacetin breath test for estimating ascorbic acid requirements. *Int J Vitam Nutr Res.* 1987;57(4):455.
9. Krumbiegel P, Günther K, Faust H, Möbius G, Hirschberg K, Schneider G. Nuclear medicine liver function tests for pregnant women and children. *Eur J Nucl Med* 1985;10:129-133.

3. Mitochondrial Liver Funktion

¹³C-Ketoisocaproic Acid Breath Test (¹³C-KICA)

■ Recommended literature

1. Parra D, Gonzalez A, Mugueta C, Martinez A, Monreal I. Laboratory approach to mitochondrial diseases. J Physiol Biochem. 2001 Sep;57(3):267-84. Review.
2. Pendtsen P, Hannestad U, Pahlsson P. Evaluation of the carbon ¹³-labeled Ketoisocaproate breath test to assess mitochondrial dysfunction in patients with high alcohol consumption. Alcohol Clin Exp Res. 1998 Nov;22(8):1792-5.
3. Lauterburg BH, Grattagliano I, Gmur R, Stalder M, Hildebrand P. Noninvasive assessment of the effect of xenobiotics on mitochondrial function in human beings: studies with acetylsalicylic acid and ethanol with the use of the carbon ¹³-labeled ketoisocaproate breath test. J Lab Clin Med. 1995 Mar;125(3):378-83.
4. Mion F, Rousseau M, Brazier JL, Minaire Y. Human hepatic macrovesicular steatosis: a noninvasive study of mitochondrial ketoisocaproic acid decarboxylation. Metabolism. 1995 Jun;44(6):699-700
5. Witschi A, Mossi S, Meyer B, Junker E, Lauterburg BH. Mitochondrial function reflected by the decarboxylation of [¹³C]ketoisocaproate is impaired in alcoholics. Alcohol Clin Exp Res. 1994 Aug;18(4):951-5.

4. Cytosolic Liver Funktion ¹³C-Galactose Breath Test

■ Recommended literature

1. Herold C, Regn S, Ganslmayer M, Ocker M, Hahn EG, Schuppan D. Can quantitative tests of liver function discriminate between different etiologies of liver cirrhosis? *Dig Dis Sci.* 2002 Dec;47(12):2669-73.
2. Herold C, Ganslmayer M, Deynet C, Hahn EG, Schuppan D. Quantitative testing of liver function compared to prognostic scores in patients with primary biliary cirrhosis. *Liver.* 2002 Apr;22(2):159-65.
3. Suzuki S, Ishii Y, Asai S, Kohno T, Mazaki T, Takahashi Y, Iwai S, Ishikawa K. 1-(¹³C) breath test of galactose and fructose for quantitative liver function. *J Surg Res.* 2001 Mar;96(1):90-5.
4. Suzuki S, Ishii Y, Asai S, Kohno T, Mazaki T, Takahashi Y, Iwai S, Ishikawa K. 1-(¹³C) Galactose breath test for quantitative measurement of liver function in a short period. *Digestion.* 2000;62(2-3):194-9.
5. Mion F, Geloën A, Minaire Y. Effects of ethanol and diabetes on galactose oxidative metabolism and elimination in rats. *Can J Physiol Pharmacol.* 1999 Mar;77(3):182-7.
6. Mion F, Rousseau M, Scoazec JY, Berger F, Minaire Y. [¹³C]-Galactose breath test: correlation with liver fibrosis in chronic hepatitis C. *Eur J Clin Invest.* 1999 Jul;29(7):624-9.
7. Becker M. ¹³C breath test for measurement of liver function. *Gut.* 1998 Nov;43 Suppl 3:S25-7.
8. Shreeve WW, Shoop JD, Ott DG, McInteer BB. Test for alcoholic cirrhosis by conversion of [¹⁴C]- or [¹³C]galactose to expired CO₂. *Gastroenterology.* 1976 Jul;71(1):98-101.