


# Protocols - $^{13}\text{C}$ Breath Tests - PANCREAS



## Table of Contents

Introduction.....	3
• <sup>13</sup> C Breath Test protocols .....	3
• <sup>13</sup> C Breath Testing: principle and requirements.....	3
• Preparation of the patient.....	3
• Administration of <sup>13</sup> C labeled substrate .....	3
• Collection of breath samples .....	3
• Measurement of <sup>13</sup> C enrichment.....	3
• Calculation of the end result.....	3
• Applications .....	4
• Literature .....	4
• Note.....	4
• Authors.....	4
Test Protocols Pancreas .....	5
• Lipase Activity .....	5
<sup>13</sup> C-Mixed Triglyceride Breath Test	
• Exocrine Pancreatic Insufficiency.....	6
<sup>13</sup> C-Cholesteryl Octanoate Breath Test	
Literature .....	7
Pancreas .....	7
• Lipase Activity .....	7
<sup>13</sup> C-Mixed Triglyceride Breath Test	
• Exocrine Pancreatic Insufficiency.....	9
<sup>13</sup> C-Cholesteryl Octanoate Breath Test	

## Introduction

### ■ **<sup>13</sup>C Breath Test protocols**

This folder contains a set of protocols describing the principles and general test procedures for today's most relevant <sup>13</sup>C Breath Tests to study specific functions of the liver, pancreas, stomach and intestine. The list will be updated regularly adding additional tests or additional information on the already described tests.

The information is meant as a start to enter the field of stable isotope <sup>13</sup>C Breath testing initiated by the interest in a specific test. To actually be able to introduce a test in your hospital you must familiarize yourself with basic knowledge of breath testing with <sup>13</sup>C substrates and the existing knowledge on the particular application of interest. There is no such thing as a standard protocol for all tests.

### ■ **<sup>13</sup>C Breath Testing: principle and requirements**

A <sup>13</sup>C Breath Test consists of the administration to a patient of a <sup>13</sup>C labeled substrate that is metabolized by a specific enzyme system resulting in <sup>13</sup>CO<sub>2</sub> as the end product. To monitor the enzyme response <sup>13</sup>C enrichment is measured in breath CO<sub>2</sub>.

The total procedure of <sup>13</sup>C Breath testing includes the definition of the preparation of the patient before the test, administration of the <sup>13</sup>C labeled substrate, collection of breath samples, measurement of <sup>13</sup>C enrichment in breath CO<sub>2</sub> and the calculation of an end result.

### ■ **Preparation of the patient**

In general, tests will be performed in the fasting state and the patient should be at a low and stable level of natural <sup>13</sup>C abundance. Therefore, the patient must be instructed to avoid eating <sup>13</sup>C enriched foods such as corn products, cane sugar, pineapple and tequila the last days before the test and to come to the clinic fasted. In certain cases (<sup>13</sup>C Lactose-Ureide breath test) the patient must be pretreated with unlabeled substrate to stimulate the involved enzyme system.

### ■ **Administration of <sup>13</sup>C labeled substrate**

The test substrate may be administered as a simple solution in water with or without a standardized test meal. Sometimes it needs to be incorporated into a specific ingredient of the meal. The test meal and the dose of substrate may be different for adults and children.

### ■ **Collection of breath samples**

Every protocol has its own time schedule of breath collections. The number of samples may be as small as 2 or more than 20. To define the <sup>13</sup>C enrichment in breath CO<sub>2</sub>, it is also necessary to obtain at least two breath samples before the ingestion of the <sup>13</sup>C substrate to determine the natural background of <sup>13</sup>C abundance. The methodology of collecting breath samples is dependent on the technology to determine the <sup>13</sup>C enrichment. The protocols are based on Continuous Flow Isotope Ratio Mass Spectrometry as the analytical technique. In this case breath samples are simply blown through a straw into special 10 ml gas collection tubes that directly fit into the sample tray of the instruments. In case of Infrared technology special bags provided by the instrument manufacturer must be used.

### ■ **Measurement of <sup>13</sup>C enrichment**

To determine the <sup>13</sup>C abundance in breath CO<sub>2</sub> you need the availability of Isotope Ratio Mass Spectrometry (IRMS) or specialized Infrared instrumentation. The protocols are based on Isotope Ratio Mass Spectrometry. For a number of tests (Aminopyrine, Methacetin, Urea) Infrared Spectroscopy has proven to be a valid alternative analytical technique. For other tests Infrared technology has not yet been validated so far. In principle the test substrate is not a determinant of the validity of the analytical technique. It is the level of <sup>13</sup>C enrichment that determines the analytical requirement. Validation of Infrared analysis for other application is recommended, as it is recommended to validate any breath test in your own clinical laboratory. You may have instrumentation available or contact a service center for the analyses.

### ■ **Calculation of the end result**

For some tests the only calculation needed is the subtraction of the natural background value from the measured value at a defined time. In other cases it is necessary to calculate the amount of <sup>13</sup>C that is recovered in breath during the experimental period. In a third type of application the time course of the enrichment appearance is of importance requiring calculation of the rate of appearance.

## ■ Applications

In the present update the following tests have been described:

### Table of contents of applications

Function	<sup>13</sup> C Substrate	Page
<b>Pancreas</b>		
1. Lipase Activity	<sup>13</sup> C-Mixed Triglyceride	5
2. Exocrine Pancreatic Insufficiency	<sup>13</sup> C-Cholesteryl-Octanoate	6
<b>Literature</b>		7

## ■ Literature

Included is a list of literature references that will introduce you to the most important articles describing aspects of the different tests described in the protocols.

## ■ Note

Great care has been taken over the composition of the text, figures and tables. The possibility of errors however, cannot be excluded completely. Therefore Campro Scientific GmbH and the authors cannot accept any legal or other liability with respect to incorrect details and their consequences. The authors would be grateful to receive suggestions for improvements and information about errors.

This information may not be photocopied, duplicated or translated into another language neither partially nor completely, in any form without the written agreement of Campro Scientific GmbH.

## ■ Authors

1. Dr. F. Stellaard  
University Hospital Groningen, the Netherlands  
Dept. Pathology and Laboratory Medicine  
And Centre for Liver, Intestinal and Metabolic Disease

2. Dr. Ahmad Rajabi  
Campro Scientific GmbH  
Berlin, Germany

© Campro Scientific GmbH

### European Headquarters

P.O. Box 37 03 31  
D-14133 Berlin  
Germany  
Tel. : +49(0)30.629.01.89.0  
Fax : +49(0)30.629.01.89.89  
E-mail : [info@campro.eu](mailto:info@campro.eu)  
Web : [www.campro.eu](http://www.campro.eu)

### Dutch Sales Office

P.O. Box 316  
NL-3900 AH Veenendaal  
The Netherlands  
Tel. : +31(0)318.529.437  
Fax : +31(0)318.542.181  
E-mail : [info@campro.eu](mailto:info@campro.eu)  
Web : [www.campro.eu](http://www.campro.eu)

*Please contact us for technical and price information.*

## Test Protocols Pancreas

### 1. Lipase Activity

#### <sup>13</sup>C-Mixed Triglyceride Breath Test

##### ■ Principle

<sup>13</sup>C-Mixed Triglyceride consists of a Triglyceride containing 2 Stearic Acid molecules at the sn-1 and sn-3 position and one Octanoic Acid molecule at the sn-2 position. The Octanoic Acid molecule is labeled with <sup>13</sup>C at the carboxyl carbon. The official terminology is 1,3-Distearyl-2-<sup>13</sup>C-Carboxy-Octanoyl Glycerol. After administration of <sup>13</sup>C-Mixed Triglyceride with a meal, <sup>13</sup>C -Mixed Triglyceride will be transported to the small intestine where it is digested by pancreatic lipase. After digestion <sup>13</sup>C-Octanoic Acid Monoglyceride is effectively absorbed and oxidized, preferentially in the liver. The rate of appearance of <sup>13</sup>C in breath CO<sub>2</sub> is considered a reflection of the rate of digestion, i.e. pancreatic lipase activity.

##### ■ Applicability of <sup>13</sup>C-Mixed Triglyceride Breath Test

<sup>13</sup>C-Mixed Triglyceride Breath Test has so far been applied to adults, children and neonates.

##### ■ Applications

<sup>13</sup>C-Mixed Triglyceride Breath Test is used to detect diminished pancreatic lipase activity in the small intestine by comparison with a control range. The result reflects overall lipase activity. The response of patients to pancreatic lipase supplementation can be documented using the <sup>13</sup>C-Mixed Triglyceride Breath Test.

##### ■ Protocol

**Adults:** The <sup>13</sup>C -Mixed Triglyceride Breath Test is performed after an overnight fast. A dose of 5 mg/kg body weight of <sup>13</sup>C-Mixed Triglyceride is administered. For this purpose the <sup>13</sup>C-Mixed Triglyceride is mixed into habitual breakfast or preferentially into a standardized breakfast. The meal may consist of 100 g toast with 0.25 g butter per kg body weight. The <sup>13</sup>C-Mixed Triglyceride is mixed into the butter. Breath samples are collected before (2x) and every 30 minutes for 360 minutes (6 h) after ingestion of the <sup>13</sup>C-Mixed Triglyceride. <sup>13</sup>C enrichment in breath CO<sub>2</sub> is determined by Isotope Ratio Mass Spectrometry (IRMS). The cumulative percentage of <sup>13</sup>C recovered in breath during the 360 minute collection period is used as the diagnostic parameter.

**Children:** A dose of 10-20 mg/kg body weight is used for children below 30 kg body weight, 5mg/kg is applied to children with body weight of 30 kg over. The dose is administered in a liquid test meal (Calogen, Nutri-mix) giving about 0.7 g/kg of fat. For neonates the label is administered directly into the mouth dissolved in olive oil during the milk feeding. The same time schedule for breath collections can be used. The sample collection technique must be adapted to the age of the children.

##### ■ Interpretation of test results

It is advised to obtain your own internal control values. Generally a cut-off value of 22% for the 6 hours cumulative recovery may serve as a starting point in the case of adults and children. Values below 22% indicate for a diminished lipase activity.

##### ■ Precautions

No side effects have been reported so far. It needs to be mentioned that in exceptional cases the absorption of long chain fatty acids may be disturbed leading to fat malabsorption. This will not affect the results of the <sup>13</sup>C-Mixed Triglyceride Breath Test.

##### ■ Summary

	Dose	Samples	
Adults	5 mg/kg body weight <sup>13</sup> C -Mixed Triglyceride	2	Before administration
		12	Every 30 minutes for 360 minutes after administration (6 hours)
Children	10-20 mg/kg body weight for children below 30 kg body weight, 5 mg/kg body weight for children above 30 kg body weight		The sample collection technique has to be adapted to the age of the child

## 2. Exocrine Pancreatic Insufficiency <sup>13</sup>C-Cholesteryl Octanoate Breath Test

### ■ Principle

(1-<sup>13</sup>C)-Cholesteryl Octanoate is a natural cholesterol ester labeled with the non-radioactive isotope <sup>13</sup>C in the Octanoic Acid part of the molecule. After oral administration Cholesteryl (1-<sup>13</sup>C) Octanoate is hydrolyzed by pancreatic cholesterol esterase. <sup>13</sup>C-Octanoic Acid is released, absorbed and effectively oxidized primarily in the liver. The degree of appearance of <sup>13</sup>C in breath CO<sub>2</sub> reflects the degree of pancreatic cholesterol esterase activity and therewith pancreatic function.

### ■ Applicability of <sup>13</sup>C-Cholesteryl Octanoate Breath Test

<sup>13</sup>C-Cholesteryl Octanoate Breath Test has so far been applied to adults.

### ■ Applications

<sup>13</sup>C-Cholesteryl Octanoate Breath Test is used to detect exocrine pancreatic insufficiency by comparison with a control range.

### ■ Protocol

**Adults:** The <sup>13</sup>C-Cholesteryl Octanoate Breath Test is performed after an overnight fast. A dose of 500 mg (1-<sup>13</sup>C)-Cholesteryl Octanoate together with 800 mg of unlabeled material is administered orally after dissolution in a liquid meal containing emulsifiers such as glycerol, lecithin and vegetable broth in order to solubilize the hydrophobic substrate. Breath samples are collected before (2x) and every 15 minutes for 180 minutes (3 h) after ingestion of the (1-<sup>13</sup>C)-Cholesteryl Octanoate. <sup>13</sup>C enrichment in breath CO<sub>2</sub> is determined by Isotope Ratio Mass Spectrometry (IRMS). The cumulative percentage of <sup>13</sup>C recovered in breath during the 180 minutes collection period is used as the diagnostic parameter

### ■ Interpretation of test results

It is advised to obtain your own internal control values. Generally a cut-off value of 6% for the 3 hours of cumulative recovery may serve as a starting point in the case of adults. Values below 6% indicate for a diminished pancreatic function.

### ■ Precautions

No precautions have been described so far.

### ■ Summary

	Dose	Samples	
		500 mg Cholesteryl (1- <sup>13</sup> C) Octanoate + 800 mg unlabeled material	2
		12	Every 15 minutes for 180 minutes after administration (3 h)



## Literature

### Pancreas

#### 1. Lipase Activity

##### <sup>13</sup>C-Mixed Triglyceride Breath Test

#### ■ Recommended literature

1. van Dijk-van Aalst K, Van Den Driessche M, van Der Schoor S, Schiffelers S, van't Westeinde T, Ghoos Y, Veereman-Wauters G. <sup>13</sup>C mixed triglyceride breath test: a noninvasive method to assess lipase activity in children. *J Pediatr Gastroenterol Nutr.* 2001 May;32(5):579-85.
2. Ling SC, Amarri S, Slater C, Hollman AS, Preston T, Weaver LT. Liver disease does not affect lipolysis as measured with the <sup>13</sup>C-mixed triacylglycerol breath test in children with cystic fibrosis. *J Pediatr Gastroenterol Nutr.* 2000 Apr;30(4):368-72.
3. Boedeker C, Goetze O, Pfaffenbach B, Luypaerts A, Geypens B, Adamek RJ. <sup>13</sup>C mixed-triglyceride breath test: isotope selective non-dispersive infrared spectrometry in comparison with isotope ratio mass spectrometry in volunteers and patients with chronic pancreatitis. *Scand J Gastroenterol.* 1999 Nov;34(11):1153-6.
4. Kalivianakis M, Elstrodt J, Havinga R, Kuipers F, Stellaard F, Sauer PJ, Vonk RJ, Verkade HJ. Validation in an animal model of the carbon <sup>13</sup>-labeled mixed triglyceride breath test for the detection of intestinal fat malabsorption. *J Pediatr.* 1999 Oct;135(4):444-50.
5. Wutzke KD, Radke M, Breuel K, Gurk S, Lafrenz JD, Heine WE. Triglyceride oxidation in cystic fibrosis: a comparison between different <sup>13</sup>C-labeled tracer substances. *J Pediatr Gastroenterol Nutr.* 1999 Aug;29(2):148-54.
6. Manson WG, Coward WA, Harding M, Weaver LT. Development of fat digestion in infancy. *Arch Dis Child Fetal Neonatal Ed.* 1999 May;80(3):F183-7.
7. Perri F, Pastore M, Festa V, Clemente R, Quitadamo M, D'Altilia MR, Niro G, Paolucci P, Andriulli A. Intraduodenal lipase activity in celiac disease assessed by means of <sup>13</sup>C mixed-triglyceride breath test. *J Pediatr Gastroenterol Nutr.* 1998 Oct;27(4):407-10.
8. De Boeck K, Delbeke I, Eggermont E, Veereman-Wauters G, Ghoos Y. Lipid digestion in cystic fibrosis: comparison of conventional and high-lipase enzyme therapy using the mixed-triglyceride breath test. *J Pediatr Gastroenterol Nutr.* 1998 Apr;26(4):408-11.
9. Loser C, Brauer C, Aygen S, Hennemann O, Folsch UR. Comparative clinical evaluation of the <sup>13</sup>C-mixed triglyceride breath test as an indirect pancreatic function test. *Scand J Gastroenterol.* 1998 Mar;33(3):327-34.
10. Swart GR, Baartman EA, Wattimena JL, Rietveld T, Overbeek SE, van den Berg JW. Evaluation studies of the <sup>13</sup>C-mixed triglyceride breath test in healthy controls and adult cystic fibrosis patients with exocrine pancreatic insufficiency. *Digestion.* 1997;58(5):415-20.
11. Kalivianakis M, Verkade HJ, Stellaard F, van der Were M, Elzinga H, Vonk RJ. The <sup>13</sup>C-mixed triglyceride breath test in healthy adults: determinants of the <sup>13</sup>CO<sub>2</sub> response. *Eur J Clin Invest.* 1997 May;27(5):434-42.
12. Amarri S, Harding M, Coward WA, Evans TJ, Weaver LT. <sup>13</sup>Carbon mixed triglyceride breath test and pancreatic enzyme supplementation in cystic fibrosis. *Arch Dis Child.* 1997 Apr;76(4):349-51.
13. Maes BD, Ghoos YF, Geypens BJ, Hiele MI, Rutgeerts PJ. Relation between gastric emptying rate and rate of intraluminal lipolysis. *Gut.* 1996 Jan;38(1):23-7.

14. Vantrappen GR, Rutgeerts PJ, Ghoo YF, Hiele MI. Mixed triglyceride breath test: a noninvasive test of pancreatic lipase activity in the duodenum. *Gastroenterology*. 1989 Apr;96(4):1126-34.
15. Pomposelli JJ, Moldawer LL, Palombo JD, Babayan VK, Bistran BR, Blackburn GL. Short-term administration of parenteral glucose-lipid mixtures improves protein kinetics in portacaval shunted rats. *Gastroenterology*. 1986 Aug;91(2):305-12.
16. Ghoo YF, Vantrappen GR, Rutgeerts PJ, Schurmans PC. A mixed-triglyceride breath test for intraluminal fat digestive activity. *Digestion*. 1981;22(5):239-47.



## 2. Exocrine Pancreatic Insufficiency <sup>13</sup>C-Cholesteryl Octanoate Breath Test

### ■ Recommended literature

1. Ventrucchi M, Cipolla A, Ubalducci GM, Roda A, Roda E. <sup>13</sup>C labelled Cholesteryl Octanoate breath test for assessing pancreatic exocrine insufficiency. Gut. 1998 Jan;42(1):81-7.
2. Bruno MJ, Hoek FJ, Delzenne B, van Leeuwen DJ, Schteingart CD, Hofmann AF, Tytgat GN. Simultaneous assessments of exocrine pancreatic function by cholesteryl-[<sup>14</sup>C]octanoate breath test and measurement of plasma p-aminobenzoic acid. Clin Chem. 1995 Apr;41(4):599-604.
3. Mundlos S, Kuhnelt P, Adler G. Monitoring enzyme replacement treatment in exocrine pancreatic insufficiency using the cholesteryl octanoate breath test. Gut. 1990 Nov;31(11):1324-8.
4. Cole SG, Rossi S, Stern A, Hofmann AF. Cholesteryl octanoate breath test. Preliminary studies on a new noninvasive test of human pancreatic exocrine function. Gastroenterology. 1987 Dec;93(6):1372-80.
5. Mundlos S, Rhodes JB, Hofmann AF. The Cholesteryl Octanoate breath test: a new procedure for detection of pancreatic insufficiency in the rat. Pediatr Res. 1987 Sep;22(3):257-61.