

D-Glucose-¹³C₆

Campro Scientific GmbH, recognized as distributor for the world's premier producer of labelled D-Glucose, is offering the largest production capacity of D-Glucose-C-d₇, D-Glucose-¹³C₆, C-d₇, and D-Glucose-¹³C₆.

Your Advatage


We continue to be the trusted source for stable isotopes for biomolecular NMR research.

Users may obtain better protein expression by using D-Glucose-¹³C₆ and Ammonium-¹⁵N Salts in combination with other biomolecular NMR Products, including:

ISOGRO™ Growth Media -¹³C; -¹³C, ¹⁵N; -¹⁵N

Amino Acid Mix -¹³C; -¹³C, ¹⁵N; -¹⁵N


Amino Acids-¹³C; -¹³C, ¹⁵N; -¹⁵N



High Resolution Solution NMR Structure of the Z. Domain of Staphylococcal Protein A. (Courtesy of M. Tashiro, R. Tejero, D.E. Zimmermann, B. Celda, B. Nilsson, G.T. Montelione, Rutgers University).

Customizing your growth media

with our D-Glucose-¹³C₆, Ammonium-¹⁵N Salts and labelled ISOGRO™ Powder or Amino Acids can be a cost-effective method to boost your expression in bacterial and yeast cultures. We can offer customized product packages to meet your formulation needs.



Quality & Reliability

you can rely on consistent, high quality products offered in gram to kilogram quantities with the fastest delivery time. Our D-Glucose-¹³C₆ is regarded as the best in the industry, with the highest isotopic enrichment and chemical purity.

Why take a risk on lower grade products?


Available from Bulk Stock

D-Glucose-¹³C₆, C-d₇

for Biomolecular NMR Studies

Campro Scientific GmbH is distributor of the world's premier producer of labelled D-Glucose, offering:

D-Glucose-C-d₇,
D-Glucose-¹³C₆, and
D-Glucose-¹³C₆,C-d₇



Specifications

Chemical Purity	Isotopic Enrichment
min. 99 %	min. 99 atom % ¹³ C 97-99 atom % D
min 99 %	min 99 atom % ¹³ C 80 atom % D
The ISOTEC Advantage	Quality Control
<p>Isotec continues to be a trusted source for reliable stable isotope labelled Biomolecular NMR products, including:</p> <p>D-Glucose-¹³C₆, ¹³C₆,C-d₇; C-d₇ Ammonium-¹⁵N Salts Growth Media-¹³C; ¹³C, ¹⁵N; ¹⁵N; ¹³C, ¹⁵N, D Algal Amino Acid Mix-¹³C; ¹³C, ¹⁵N; ¹⁵N Amino Acids-ul-¹³C, ¹⁵N Selectively Labelled Amino Acids</p> <p>With the entire production process under our control, you can count on consistent, high quality products from gram to kilogram quantities with the fastest delivery time.</p>	<p>Isotec uses controls throughout the production process and tests the final product rigorously to minimize variation and assure superior quality. Our D-Glucose-¹³C₆, C-d₇ has the highest isotopic enrichment and chemical purity available in the industry.</p>