

**NEW**  
AVAILABLE NOW!



**FANhp**

New, compact C13-urea breath test device for all suitable H. pylori test kits  
Economically priced, no PC required  
Easy to handle via touchscreen  
Connection of receipt printer possible  
Simple transfer of results via USB stick

**Features**

- Adjustable pos/neg threshold for all H. pylori breath test kits with bags
- Automatic correction of inverted samples
- Instant printing upon availability of result
- Integrated self-test
- Easy-to-read AMOLED display
- Touchscreen operable with gloves
- Convertible port connectors for various breath bags
- Service interface for software updates and maintenance
- Optional scanner for barcode operations



Europa fördert Sachsen.  
**EFRE**  
Europäischer Fonds für regionale Entwicklung



Management System  
ISO 9001:2008

www.tuv.com  
ID 9108622470

Campro Scientific GmbH - Germany

P.O. Box 37 03 31 • 14133 Berlin • Tel. +49.(0)30.629.01.89.0 • Fax +49.(0)30.629.01.89.89 • info@campro.eu • www.campro.eu

Campro Scientific GmbH - The Netherlands

P.O. Box 316 • NL-3900 AH Veenendaal • Tel. +31.(0)318.529.437 • Fax +31.(0)318.542.181 • info.nl@campro.eu • www.campro.eu