

Total Solution Sample Prep

- Extraction
- Cleanup
- Concentration

TRP System
Total
Rapid
Prep



from sample to vial



SPE
Solid
Phase
Extraction



PLE®
Pressurized
Liquid
Extraction



PowerPrep™
Multi-Column
Cleanup-System



SuperVap™
Concentration
System



CAMPRO
SCIENTIFIC

FMS
Fluid Management Systems

Total Solution Sample Prep

■ Pressurized Liquid Extraction

Solid and Semi-Solid Samples

The PLE™ is a high speed Pressurized Liquid Extraction System, designed to perform extraction of multiple samples simultaneously in minutes rather than hours, producing high recoveries and excellent precision for all analytes.

- Fish
- Feed
- Meat
- Fly Ash
- Soil
- Sludge



PLE®
Pressurized
Liquid
Extraction

■ PowerPrep™ SPE Solid Phase Extraction

Liquid Samples

The PowerPrep SPE automated One Step Extraction and Concentration System is built by design to simply, improve and increase your productivity by automating the manual steps in your sample preparation process. The PowerPrep SPE system replaces older techniques, Liquid-Liquid Extraction and outdated semi-automated instruments and automates existing manual SPE techniques.

- Drinking water
- Waste water
- Milk
- Serum
- Pharmaceutical
- Personal Care By Products



SPE
Solid
Phase
Extraction

■ PowerPrep™ GPC Sample Cleanup

The GPC Sample Cleanup System automates pesticide sample cleanup procedures specified by US EPA and FDA methods. The modularity and expansion functionality of the GPC-Prep allow the user to expand the system for use in other applications such as the separation of Dioxins and PCBs from environmental and biological samples.



GPC
Sample
Cleanup

■ PowerPrep™ Multicolumn Sample Cleanup

The PowerPrep Multi-Column Sample Cleanup System performs the sample cleanup process in less than two hours instead of days for the trace analysis of Persistent Organic Pollutants in all sample matrices.



PowerPrep™
System

■ SuperVap™ Concentration System

Direct-to-vial concentration, Standalone or integrated, Dry Bath, HEPA/Carbon Filtration for BFRs, Dioxins, Furans, PBDFs, PAHs, Pesticides and PPCPs



SuperVap™
Concentration
System

CAMPRO
SCIENTIFIC

■ Germany

Campro Scientific GmbH
Tel. +49.(0)30.629.01.89.0
info@campro.eu
www.campro.eu

■ The Netherlands

Campro Scientific GmbH
Tel. +31.(0)318.529.437
info.nl@campro.eu
www.campro.eu

FMS
Fluid Management Systems

■ USA

FMS - Fluid Management Systems, Inc.
Tel. +1.617.393.2396
sales@fms-inc.com
www.fmsenvironmental.com