



Underground Storage Tank

Our Underground Storage Tank (UST) products in water and soil matrices are purposefully designed to meet accreditation requirements for Petroleum Hydrocarbons analysis in various jurisdictions.



UST in Water PT Scheme Schedule

2026 Schedule

	Scheme #	Opens	Closes
Q	WP 372	Jan 20	Mar 6
Q	WP 375	Apr 20	Jun 4
Q	WP 378	Jul 20	Sep 3
Q	WP 381	Oct 16	Nov 30

2027 Schedule

	Scheme #	Opens	Closes
Q	WP 384	Jan 19	Mar 5
Q	WP 387	Apr 19	Jun 3
Q	WP 390	Jul 19	Sep 2
Q	WP 393	Oct 15	Nov 29



Soil (including UST in Soil) PT Schedule

2026 Schedule

	Scheme #	Opens	Closes
Q	SOIL 133	Jan 26	Mar 12
Q	SOIL 134	Apr 27	Jun 11
Q	SOIL 135	Jul 27	Sep 10
Q	SOIL 136	Oct 23	Dec 7

2027 Schedule

	Scheme #	Opens	Closes
Q	SOIL 137	Jan 25	Mar 11
Q	SOIL 138	Apr 26	Jun 10
Q	SOIL 139	Jul 26	Sep 9
Q	SOIL 140	Oct 22	Dec 6

Schedule subject to change - see Waters ERA's website at eraqc.com

For the latest products and information, please visit us online at eraqc.com

Contents

CRM Certified Reference Material
PT Proficiency Testing
QR QuiK Response
RM Reference Material

All Waters ERA UST PTs open quarterly (**Q**) or biannually (**B**), unless otherwise noted. Quarterly months are January, April, July, and October. Biannual months are January and July.

B Waters ERA NJ EPH in Soil PT opens in April and October.

Description	CRM	PT	QR	Page
Alaska BTEX in Soil	636	—	470QR	49
Alaska BTEX in Water	646	—	474QR	49
Alaska DRO in Soil	637	—	471QR	49
Alaska DRO in Water	647	—	475QR	49
Alaska GRO in Soil	635	—	469QR	49
Alaska GRO in Water	645	—	473QR	49
Alaska RRO in Soil	638	—	472QR	49
Arizona TPH in Soil	798	488 Q	798QR	49
BTEX & MTBE in Soil	761	633 Q	761QR	48
BTEX & MTBE in Water	760	643 Q	760QR	48
Diesel Range Organics in Soil	765	631 Q	765QR	48
Diesel Range Organics in Water	764	641 Q	764QR	48
Gasoline Range Organics in Soil	763	630 Q	763QR	48
Gasoline Range Organics in Water	762	640 Q	762QR	48
Massachusetts EPH in Soil	569	484 Q	569QR	50
Massachusetts VPH in Soil	568	483 Q	568QR	50
Massachusetts EPH in Water	567	482 Q	567QR	50
Massachusetts VPH in Water	566	481 B	566QR	50

Description	CRM	PT	QR	Page
New Jersey EPH in Soil	564	464 B	564QR	50
Texas High-Level Fuels in Soil	797	479 Q	797QR	49
Texas High-Level Fuels in Water	795	477 B	795QR	49
Texas Low-Level Fuels in Soil	796	478 Q	796QR	49
Texas Low-Level Fuels in Water	794	476 B	794QR	49
Total Petroleum Hydrocarbons (TPH) in Soil #1	570	632 Q	572QR	48
Total Petroleum Hydrocarbons (TPH) in Soil #2	571	632 Q	572QR	48
Total Petroleum Hydrocarbons (TPH) in Water #1	600	642 Q	602QR	48
Total Petroleum Hydrocarbons (TPH) in Water #2	601	642 Q	602QR	48
Washington HEM/SGT-HEM	519	489 Q	519QR	50
Wisconsin Gasoline Range Organics (GRO/PVOC) in Water	773	649 B	773QR	50
Wisconsin Diesel Range Organics (DRO) in Water	772	648 B	772QR	50

CRM: A reference material characterized by a metrologically valid procedure for one or more specified properties, accompanied by a reference material certificate that provides the value of the specified property, its associated uncertainty, and a statement of metrological traceability.

A complete listing of ERA's CRMs can be found on our Scope of Accreditation for general requirements for competence of reference material producers available at eraqc.com/Accreditations.

PT: A Proficiency Test (PT) is an analysis of what is often referred to as a blind sample or a sample with unknown concentrations of analytes for the purpose of evaluating a laboratory's analytical performance.

QR: Similar to a Proficiency Test, a QuiK Response (QR) is a sample with unknown concentrations. However, unlike a scheduled PT, QR is on-demand and available at any time. Plus, your results are returned within two business days. QuiK Response can be used as a bilateral PT as referenced in the IUPAC/CITAC guide: Selection and use of PT schemes for a limited number of participants – chemical analytical labs.

RM: A material, sufficiently homogeneous and stable with respect to one or more specified properties, which has been established to be fit for its intended use in a measurement process.

UST in Soil

BTEX & MTBE in Soil			
CRM	PT	Q	QR
Cat. #761	Cat. #633		Cat. #761QR

One 2 mL flame-sealed ampule requires spiking onto the ten grams of provided certified clean soil. Includes all the BTEX compounds and MTBE at 20–200 µg/kg (40–400 µg/kg for total xylenes). Use with EPA Method 8021, or other applicable methods.

Gasoline Range Organics (GRO) in Soil			
CRM	PT	Q	QR
Cat. #763	Cat. #630		Cat. #763QR

One flame-sealed ampule with 20 g of soil spiked with unleaded regular gasoline in the range 100–2000 mg/kg. Use with purge and trap and modified EPA Method 8015, or other applicable GC/FID methods. Also use to test for BTEX in gasoline.

Note: This standard is not compliant with the NELAC concentration ranges for the BTEX analytes. If a NELAC-compliant sample for these analytes is required, use Volatiles in Soil, Cat. #623 or BTEX & MTBE Soil, Cat. #633.

Diesel Range Organics (DRO) in Soil			
CRM	PT	Q	QR
Cat. #765	Cat. #631		Cat. #765QR

One flame-sealed ampule with 20 g of soil spiked with #2 Diesel Fuel in the range 300–3000 mg/kg. Use with modified EPA Method 8015, or other applicable GC/FID methods.

Total Petroleum Hydrocarbons (TPH) in Soil #1			
CRM	PT	Q	QR
Cat. #570	Cat. #632		Cat. #572QR

One screw-top bottle with 50 g of soil to be analyzed for total petroleum hydrocarbons (TPH). Use with EPA IR, Gravimetric Methods 8440 and 9071B.

Non-polar extractable material (TPH) (Gravimetric).....300–3000 mg/kg
 Non-polar extractable material (TPH) (IR).....300–3000 mg/kg

Total Petroleum Hydrocarbons (TPH) in Soil #2			
CRM	PT	Q	QR
Cat. #571	Cat. #632		Cat. #572QR

One screw-top bottle contains 50 g of soil with TPH in the presence of interfering fatty acids. Use with EPA Methods 8440, 9071B.

Non-polar extractable material (TPH) (Gravimetric).....300–3000 mg/kg
 Non-polar extractable material (TPH) (IR).....300–3000 mg/kg



UST in Water

BTEX & MTBE in Water			
CRM	PT	Q	QR
Cat. #760	Cat. #643		Cat. #760QR

One 2 mL flame-sealed ampule yields in excess of 200 mL after dilution. Use with EPA Methods 602, 8021, or other applicable methods. Includes all BTEX compounds and MTBE at 5–300 µg/L after dilution.

Gasoline Range Organics (GRO) in Water			
CRM	PT	Q	QR
Cat. #762	Cat. #640		Cat. #762QR

One 2 mL flame-sealed ampule yields up to 2 liters after dilution. Use with both purge and trap, and modified EPA Method 8015, or other applicable GC/FID methods to test for GRO at 400–4000 µg/L. Also use to test for BTEX in gasoline.

Diesel Range Organics (DRO) in Water			
CRM	PT	Q	QR
Cat. #764	Cat. #641		Cat. #764QR

One 2 mL flame-sealed ampule yields up to 2 liters after dilution. Use with modified EPA Method 8015, or other applicable GC/FID methods. Includes #2 Diesel Fuel at 800–6000 µg/L.

Total Petroleum Hydrocarbons (TPH) in Water #1			
CRM	PT	Q	QR
Cat. #600	Cat. #642		Cat. #602QR

One liter whole-volume bottle is ready to analyze for total petroleum hydrocarbons (TPH) without interfering fatty acids. Use with EPA Methods 418.1, 1664, 5520.

Total petroleum hydrocarbons.....20–200 mg/L

Total Petroleum Hydrocarbons (TPH) in Water #2			
CRM	PT	Q	QR
Cat. #601	Cat. #642		Cat. #602QR

One liter whole-volume bottle is ready to analyze for TPH in water in the presence of interfering fatty acids. Use with EPA Methods 418.1, 1664, 5520, 8440.

Total petroleum hydrocarbons.....20–200 mg/L



Learn more about
Underground Storage products

Alaska UST in Water

Alaska GRO in Water

CRM	PT		QR
Cat. #645	—	—	Cat. #473QR

One 2 mL flame-sealed ampule. Use with method AK101 for unleaded regular gasoline at 100–500 µg/L after dilution.

Alaska DRO in Water

CRM	PT		QR
Cat. #647	—	—	Cat. #475QR

One 2 mL flame-sealed ampule. Use with method AK102 for #2 Diesel Fuel at 800–2300 µg/L after dilution.

Alaska BTEX in Water

CRM	PT		QR
Cat. #646	—	—	Cat. #474QR

One 2 mL flame-sealed ampule. Use with method AK101 for all BTEX analytes at 5–30 µg/L after dilution.

Alaska UST in Soil

Alaska GRO in Soil

CRM	PT		QR
Cat. #635	—	—	Cat. #469QR

One 20 mL flame-sealed ampule with 10 g of soil and 10 mL of methanol with unleaded regular gasoline at 30–1500 mg/kg. Use with method AK101.

Alaska DRO in Soil

CRM	PT		QR
Cat. #637	—	—	Cat. #471QR

One flame-sealed ampule with 20 g of soil spiked with #2 Diesel Fuel at 30–1500 mg/kg. Use with method AK102.

Alaska RRO in Soil

CRM	PT		QR
Cat. #638	—	—	Cat. #472QR

One flame-sealed ampule with 20 g of soil with Residual Range Organic fuels at 150–2000 mg/kg. Use with method AK103.

Alaska BTEX in Soil

CRM	PT		QR
Cat. #636	—	—	Cat. #470QR

One 2 mL flame-sealed ampule along with clean soil matrix for spiking. Use with method AK101 for all BTEX analytes at 5–100 mg/kg after spiking.

Arizona UST in Soil

Arizona TPH in Soil

CRM	PT	Q	QR
Cat. #798	Cat. #488		Cat. #798QR

One ready-to-use flame-sealed ampule with 30 g of soil with Oil Range Organics and #2 Diesel Fuel. Use with method 8015AZ for TPH in the range 300–400 mg/kg. Also includes two carbon ranges.

Texas TPH in Water

All Texas TPH PT standards are designed for use with TNRCC 1005 method. The standards meet the requirements of all states that accredit for these methods including Texas, Louisiana, and Oklahoma.

Texas Low-Level Fuels (TPH) in Water

CRM	PT	B	QR
Cat. #794	Cat. #476		Cat. #794QR

One 2 mL flame-sealed ampule yields in excess of 200 mL after dilution. Contains unleaded regular gasoline and #2 Diesel Fuel resulting in TPH in the range 5–10 mg/L.

Texas High-Level Fuels (TPH) in Water

CRM	PT	B	QR
Cat. #795	Cat. #477		Cat. #795QR

One 2 mL flame-sealed ampule yields in excess of 200 mL after dilution. Contains unleaded regular gasoline and #2 Diesel Fuel resulting in TPH in the range 20–100 mg/L.

Texas TPH in Soil

Texas Low-Level Fuels (TPH) in Soil

CRM	PT	Q	QR
Cat. #796	Cat. #478		Cat. #796QR

One ready-to-use flame-sealed ampule with 20 g of soil with unleaded gasoline and #2 Diesel Fuel for TPH in the range 50–100 mg/kg.

Texas High-Level Fuels (TPH) in Soil

CRM	PT	Q	QR
Cat. #797	Cat. #479		Cat. #797QR

One ready-to-use flame-sealed ampule with 20 g of soil with unleaded gasoline and #2 Diesel Fuel for TPH in the range 1000–20,000 mg/kg.

Wisconsin GRO/PVOC/DRO Method UST

All Wisconsin UST PT standards are designed for use with Wisconsin GRO/PVOC or DRO Methods. The standards meet the requirements of all states that accredit for these methods including Wisconsin and Minnesota.

Wisconsin Gasoline Range Organics (GRO/PVOC) in Water

CRM Cat. #773	PT Cat. #649	B	QR Cat. #773QR
-------------------------	------------------------	----------	--------------------------

One 2 mL flame-sealed ampule yields in excess of 200 mL after dilution. Includes ten gasoline range synthetic organic compounds as defined by Wisconsin. Use with Wisconsin GRO/PVOC Method.

Wisconsin Diesel Range Organics (DRO) in Water

CRM Cat. #772	PT Cat. #648	B	QR Cat. #772QR
-------------------------	------------------------	----------	--------------------------

One 2 mL flame-sealed ampule yields up to 2 liters after dilution. Includes ten diesel range synthetic organic compounds in the range 200–600 µg/L. Use with the Wisconsin DRO Method.

Washington HEM/SGT-HEM Method UST

The Washington UST PT standard is designed for use with EPA Method 1664 for HEM/SGT-HEM.

Washington HEM/SGT-HEM

CRM Cat. #519	PT Cat. #489	Q	QR Cat. #519QR
-------------------------	------------------------	----------	--------------------------

One 5 mL flame-sealed ampule yields up to 2 liters after dilution. Use with EPA Method 1664 to measure HEM/SGT-HEM at 5–100 mg/L.

New Jersey EPH

The New Jersey EPH in Soil standard is designed for use with the NJ Extractable Petroleum Hydrocarbons Method.

New Jersey EPH in Soil

CRM Cat. #564	PT Cat. #464	B	QR Cat. #564QR
-------------------------	------------------------	----------	--------------------------

One flame-sealed ampule with 20 g soil containing EPH in the range of 300–3000 mg/kg.

B The NJ EPH in Soil PT studies open in April and October.

Massachusetts Hydrocarbons in Water

All Massachusetts UST PT standards are designed for use with Massachusetts Volatile Petroleum Hydrocarbon or Extractable Petroleum Hydrocarbon Methods. The standards meet the requirements of all states that accredit for these methods including Massachusetts, North Carolina, and Washington when reporting the Massachusetts carbon ranges.

Massachusetts VPH in Water

CRM Cat. #566	PT Cat. #481	B	QR Cat. #566QR
-------------------------	------------------------	----------	--------------------------

One 2 mL flame-sealed ampule yields in excess of 200 mL after dilution. Contains volatile petroleum hydrocarbon fuels (VPH) in the range 400–4000 µg/L. Use with the Massachusetts Volatile Petroleum Hydrocarbon Method for multiple carbon ranges, BTEX compounds and MTBE.

Massachusetts EPH in Water

CRM Cat. #567	PT Cat. #482	Q	QR Cat. #567QR
-------------------------	------------------------	----------	--------------------------

One 2 mL flame-sealed ampule yields up to 2 liters after dilution. Contains extractable petroleum hydrocarbon fuels (EPH) in the range 800–6000 µg/L. Use with the Massachusetts Extractable Petroleum Hydrocarbon Method for multiple carbon ranges and PAH compounds.

Massachusetts Hydrocarbons in Soil

Massachusetts VPH in Soil

CRM Cat. #568	PT Cat. #483	Q	QR Cat. #568QR
-------------------------	------------------------	----------	--------------------------

One flame-sealed ampule with 20 g soil with VPH fuels. Contains volatile petroleum hydrocarbon fuels (VPH) in the range 100–2000 mg/kg. Use with the Massachusetts Volatile Petroleum Hydrocarbon Method for multiple carbon ranges, BTEX compounds and MTBE.

Massachusetts EPH in Soil

CRM Cat. #569	PT Cat. #484	Q	QR Cat. #569QR
-------------------------	------------------------	----------	--------------------------

One flame-sealed ampule with 20 g soil with EPH fuels. Contains extractable petroleum hydrocarbon fuels (EPH) in the range 300–3000 mg/kg. Use with the Massachusetts Extractable Petroleum Hydrocarbon Method for multiple carbon ranges and PAH compounds.

We focus on quality and service, so you can focus on your business.

Unmatched Technical Expertise

As your Partner in Quality, our goal is to help you maintain successful PT performance, solve routine analysis challenges, and improve corrective actions. Whether it's organic and/or inorganic chemistry, microbiology, analytical instrumentation or methods, our experts are ready to help you with:

- Method interpretations
- Prep and analytical questions
- Instrumentation troubleshooting
- Quality control issues
- Calibration issues

World-Class Customer Service

Our customer service team understands that you are faced with a myriad of requirements to maintain your laboratory accreditation. Each of our representatives has helped solve questions from customers with the same types of challenges. Your dedicated customer service representative has the experience and knowledge to help you through every step of the process.

For more information, contact our customer service team at 800.372.0122 / +1.303.431.8454. or email at era_info@waters.com.



CONTACT US
eraqc.com/contact

